



Anti-NRK (aa 1483-1582) polyclonal antibody (DPAB-DC862)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	The mouse ortholog of this gene encodes a protein kinase required for JNK activation. The encoded protein may be involved in the induction of actin polymerization in late embryogenesis.[provided by RefSeq, Jun 2010]
Immunogen	NRK (NP_940867, 1483 a.a. ~ 1582 a.a) partial recombinant protein with GST tag. The sequence is VEANEQLFKKILEMWKDIPSSIAFEQTQRTTGWGQKAIEVRSLSRVLESELKRRSIKKL RFLCTRGDKLFFTSTLRNHHSRVYFMTLGKLEELQSNDYDV
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name [NRK Nik related kinase \[Homo sapiens \(human\) \]](#)

Official Symbol	NRK
Synonyms	NRK; Nik related kinase; NESK; nik-related protein kinase;
Entrez Gene ID	203447
Protein Refseq	NP_940867
UniProt ID	Q7Z2Y5
Chromosome Location	Xq22.3
Pathway	TNF receptor signaling pathway;
Function	ATP binding; receptor signaling protein serine/threonine kinase activity; small GTPase regulator activity;