



Anti-ERH (full length) polyclonal antibody (DPAB-DC893)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	ERH (enhancer of rudimentary homolog (Drosophila)) is a protein-coding gene. Diseases associated with ERH include retinal perforation, and dental fluorosis.
Immunogen	ERH (AAH14301, 1 a.a. ~ 104 a.a) full-length recombinant protein with GST tag. The sequence is MSHTILLVQPTKRPEGRITYADYESVNECMGVCKMYEEHLKRMNPNSPSITYDISQLFDF IDDLADLSCLVYRADTQTYQPYNKDWIKEKIYVLLRRQAQQAGK
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	ERH enhancer of rudimentary homolog (Drosophila) [Homo sapiens (human)]
Official Symbol	ERH

Synonyms	ERH; enhancer of rudimentary homolog (Drosophila); DROER; enhancer of rudimentary homolog;
Entrez Gene ID	2079
Protein Refseq	NP_004441
UniProt ID	A0A024R6D4
Chromosome Location	14q24.1
Function	poly(A) RNA binding;