



Anti-YEATS2 polyclonal antibody (DPAB-DC898)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	YEATS2 (YEATS domain containing 2) is a protein-coding gene. GO annotations related to this gene include TBP-class protein binding and protein binding. An important paralog of this gene is YEATS4.
Specificity	Specific to recombinant protein GX0258. This antibody detects mYEATS2 protein.
Immunogen	Recombinant GST fusion protein corresponding to 211 amino acids of mouse Yeats2.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Conjugate	Unconjugated
Applications	WB,
Format	Liquid
Size	50 µg
Buffer	In PBS (50% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	Yeats2 YEATS domain containing 2 [Mus musculus (house mouse)]
Official Symbol	YEATS2
Synonyms	YEATS2; YEATS domain containing 2; BC042768; mKIAA1197; YEATS domain-containing protein 2;
Entrez Gene ID	208146
Protein Refseq	NP_001028409
UniProt ID	B9EKJ4
Chromosome Location	16 A3; 16
Pathway	Chromatin modifying enzymes; HATs acetylate histones.
Function	TBP-class protein binding;