



Anti-KLHL29 polyclonal antibody (DPAB-DC900)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	KLHL29 (kelch-like family member 29) is a protein-coding gene. An important paralog of this gene is KLHL21.
Specificity	Specific to recombinant protein GX0065. This antibody detects mKLHL29 protein.
Immunogen	Recombinant GST fusion protein corresponding to 278 amino acids of mouse Klhl29.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Conjugate	Unconjugated
Applications	WB,
Format	Liquid
Size	50 µg
Buffer	In PBS (50% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	Klhl29 kelch-like 29 [Mus musculus (house mouse)]
-----------	---

Official Symbol	KLHL29
Synonyms	KLHL29; kelch-like 29; Gm68; Kbtbd9; mKIAA1921; A230106N14Rik; kelch-like protein 29; kelch repeat and BTB (POZ) domain containing 9; kelch repeat and BTB domain-containing protein 9;
Entrez Gene ID	208439
Protein Refseq	NP_001157965
UniProt ID	Q80T74
Chromosome Location	12 A1.1; 12
Function	molecular_function;