



# Anti-ERVK2 (aa 88-197) polyclonal antibody (DPAB-DC902)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                     |                                                                                                                                                                                                                     |
|---------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Antigen Description | DISCONTINUED: This record was withdrawn by the HUGO Gene Nomenclature Committee (HGNC).                                                                                                                             |
| Immunogen           | ERVK2 (AAB63112, 88 a.a. ~ 197 a.a) partial recombinant protein with GST tag. The sequence is<br>PQGMLNSPTICQTFVAQVLQPVRDKFSDCYIIHYVDDILCAAETRDKLIDCYTFLQTEVA<br>NAGLTIASDKIQTSTPFHYLGMQVEERKIKPQKVEIRKDTLRTLNDFKNC |
| Source/Host         | Mouse                                                                                                                                                                                                               |
| Species Reactivity  | Human                                                                                                                                                                                                               |
| Conjugate           | Unconjugated                                                                                                                                                                                                        |
| Applications        | ELISA,                                                                                                                                                                                                              |
| Size                | 50 µl                                                                                                                                                                                                               |
| Buffer              | 50 % glycerol                                                                                                                                                                                                       |
| Preservative        | None                                                                                                                                                                                                                |
| Storage             | Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.                                                                                                                                            |

## GENE INFORMATION

|                 |                                                                                        |
|-----------------|----------------------------------------------------------------------------------------|
| Gene Name       | <a href="#">ERVK2 endogenous retroviral sequence K(C4). 2 [ Homo sapiens (human) ]</a> |
| Official Symbol | ERVK2                                                                                  |

|                            |                                                                     |
|----------------------------|---------------------------------------------------------------------|
| <b>Synonyms</b>            | ERVK2; endogenous retroviral sequence K(C4), 2; HERV-K; polymerase; |
| <b>Entrez Gene ID</b>      | <a href="#">2087</a>                                                |
| <b>Chromosome Location</b> | Xq28                                                                |