



Anti-VPS8 polyclonal antibody (DPAB-DC903)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	VPS8 (vacuolar protein sorting 8 homolog (<i>S. cerevisiae</i>)) is a protein-coding gene. GO annotations related to this gene include zinc ion binding.
Specificity	Specific to recombinant protein GX0354. This antibody detects mVPS8 protein. It also recognizes human VPS8 protein.
Immunogen	Recombinant GST fusion protein corresponding to 166 amino acids of mouse Vps8.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse
Conjugate	Unconjugated
Applications	WB,
Format	Liquid
Size	50 µg
Buffer	In PBS (50% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	Vps8 vacuolar protein sorting 8 homolog (<i>S. cerevisiae</i>) [<i>Mus musculus</i> (house mouse)]
------------------	--

Official Symbol	VPS8
Synonyms	VPS8; vacuolar protein sorting 8 homolog (S. cerevisiae); AI315068; AU040738; mKIAA0804; vacuolar protein sorting-associated protein 8 homolog;
Entrez Gene ID	209018
Protein Refseq	NP_001272822
UniProt ID	Q0P5W1
Chromosome Location	16 B1; 16
Function	metal ion binding; molecular_function; zinc ion binding;