



Anti-D130043K22RIK polyclonal antibody (DPAB-DC907)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to mouse D130043K22Rik. The sequence is C-QGKIKQENKPTLH
Source/Host	Goat
Species Reactivity	Mouse
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	WB (Tissue lysate), ELISA,
Format	Liquid
Concentration	0.5 mg/mL
Size	100 µg
Buffer	In Tris saline, pH 7.3 (0.5% BSA, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	D130043K22Rik RIKEN cDNA D130043K22 gene [Mus musculus (house mouse)]
Official Symbol	D130043K22RIK

Synonyms	D130043K22RIK; RIKEN cDNA D130043K22 gene; Kiaa0319; 4930451E12Rik; dyslexia-associated protein KIAA0319 homolog;
Entrez Gene ID	210108
Protein Refseq	NP_001074520
UniProt ID	Q5SZV5
Chromosome Location	13 A3.1; 13
Function	molecular_function;