



Anti-KBTBD2 polyclonal antibody (DPAB-DC909)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	KBTBD2 (kelch repeat and BTB (POZ) domain containing 2) is a protein-coding gene. An important paralog of this gene is KLHL21.
Specificity	Specific to recombinant protein GX0633. This antibody detects mKBTBD2 protein.
Immunogen	Recombinant GST fusion protein corresponding to 233 amino acids of mouse Kbtbd2.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Conjugate	Unconjugated
Applications	WB,
Format	Liquid
Size	50 µg
Buffer	In PBS (50% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	Kbtbd2 kelch repeat and BTB (POZ) domain containing 2 [Mus musculus (house mouse)]
------------------	--

Official Symbol	KBTBD2
Synonyms	KBTBD2; kelch repeat and BTB (POZ) domain containing 2; Bklhd1; BC022962; mKIAA1489; BTB and kelch domain containing 1;
Entrez Gene ID	210973
Protein Refseq	NP_666070
UniProt ID	G3X9X1
Chromosome Location	6 B3; 6
Function	molecular_function;