



Anti-CCDC85A polyclonal antibody (DPAB-DC931)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	CCDC85A (coiled-coil domain containing 85A) is a protein-coding gene. Diseases associated with CCDC85A include neurofibromatosis. An important paralog of this gene is CCDC85B.
Specificity	Specific to recombinant protein GX1954. This antibody detects mCCDC85A protein. It also recognizes human CCDC85A protein.
Immunogen	Recombinant GST fusion protein corresponding to 157 amino acids of mouse Ccdc85a.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse
Conjugate	Unconjugated
Applications	WB,
Format	Liquid
Size	50 µg
Buffer	In PBS (50% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	Ccgc85a coiled-coil domain containing 85A [Mus musculus (house mouse)]
Official Symbol	CCDC85A
Synonyms	CCDC85A; coiled-coil domain containing 85A; coiled-coil domain-containing protein 85A;
Entrez Gene ID	216613
Protein Refseq	NP_001160133
UniProt ID	Q5SP85
Chromosome Location	11 A3.3; 11
Function	molecular_function;