



Anti-LARP4B polyclonal antibody (DPAB-DC942)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	LARP4B (La ribonucleoprotein domain family, member 4B) is a protein-coding gene. Diseases associated with LARP4B include neuronitis. GO annotations related to this gene include RNA binding. An important paralog of this gene is LARP4.
Specificity	Specific to recombinant protein GX1366. This antibody detects mLARP5 protein. It also recognizes human LARP5 protein.
Immunogen	Recombinant GST fusion protein corresponding to 193 amino acids of mouse Larp4b.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse
Conjugate	Unconjugated
Applications	WB,
Format	Liquid
Size	50 µg
Buffer	In PBS (50% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	Larp4b La ribonucleoprotein domain family, member 4B [Mus musculus (house mouse)]
Official Symbol	LARP4B
Synonyms	LARP4B; La ribonucleoprotein domain family, member 4B; Larp5; AI256361; D13Wsu64e; A630096F19; la-related protein 4B; la-related protein 5; la ribonucleoprotein domain family member 5; La ribonucleoprotein domain family, member 5;
Entrez Gene ID	217980
Protein Refseq	NP_766173
UniProt ID	Q6A0A2
Chromosome Location	13 A1; 13 3.66 cM
Function	RNA binding; nucleic acid binding; poly(A) RNA binding;