



Mouse anti-Human FBN1 polyclonal antibody (DPAB-DC953)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	This gene encodes a member of the fibrillin family. The encoded protein is a large, extracellular matrix glycoprotein that serve as a structural component of 10-12 nm calcium-binding microfibrils. These microfibrils provide force bearing structural support in elastic and nonelastic connective tissue throughout the body. Mutations in this gene are associated with Marfan syndrome, isolated ectopia lentis, autosomal dominant Weill-Marchesani syndrome, MASS syndrome, and Shprintzen-Goldberg craniosynostosis syndrome.
Immunogen	FBN1 (NP_000129, 2772 a.a. ~ 2871 a.a) partial recombinant protein with GST tag. The sequence is SNKVRILELLPALTTLTNHNRYLIESGNEDGFFKINQKEGISYLFHTKKKPVAGTYSLQI SSTPLYKKKELNQLEDKYDKDYLSGELGDNLMKMIQVLLH
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	FBN1 fibrillin 1 [Homo sapiens (human)]
Official Symbol	FBN1
Synonyms	FBN1; fibrillin 1; FBN; SGS; WMS; MASS; MFS1; OCTD; SSKS; WMS2; ACMICD; ECTOL1; GPHYSD2; fibrillin-1; fibrillin 15;
Entrez Gene ID	2200
Protein Refseq	NP_000129
UniProt ID	P35555
Chromosome Location	15q21.1
Pathway	Elastic fibre formation; Integrin cell surface interactions;
Function	calcium ion binding; extracellular matrix structural constituent; protein binding;