



Anti-USP32P2 (aa 255-354) polyclonal antibody (DPAB-DC957)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	USP32P2 (ubiquitin specific peptidase 32 pseudogene 2) is a pseudogene, and is affiliated with the lncRNA class. Diseases associated with USP32P2 include smith magenis syndrome.
Immunogen	LOC220594 (NP_665808, 255 a.a. ~ 354 a.a) partial recombinant protein with GST tag. The sequence is HIKPIYNLYAISCHSGILSGGHYVTYAKNPNCWKYCYNGSICEEHPDEIDTDSAYILFY EQQRIDYAQFLPKIDGKKMADTSSMDEDFESDYEKYCVLQ
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Cell lysate), WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	USP32P2 ubiquitin specific peptidase 32 pseudogene 2 [Homo sapiens (human)]
Official Symbol	USP32P2

Synonyms	USP32P2; ubiquitin specific peptidase 32 pseudogene 2; TL132;
Entrez Gene ID	220594
Chromosome Location	17p11.2