



Anti-CTR9 (N-terminal) polyclonal antibody (DPAB-DC963)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	CTR9 (CTR9, Paf1/RNA polymerase II complex component) is a protein-coding gene. Diseases associated with CTR9 include influenza. GO annotations related to this gene include protein binding and SH2 domain binding.
Immunogen	A synthetic peptide (conjugated with KLH) corresponding to N-terminus of mouse Ctr9.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Conjugate	Unconjugated
Applications	WB (Cell lysate), ELISA,
Format	Liquid
Size	100 µl
Buffer	In serum (0.05% sodium azide)
Preservative	0.05% Sodium Azide
Storage	Store at -20°C. For long term storage store at -80°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name [Ctr9 Ctr9, Paf1/RNA polymerase II complex component, homolog \(S. cerevisiae\) \[Mus](#)

[musculus \(house mouse\) \]](#)

Official Symbol	CTR9
Synonyms	CTR9; Ctr9, Paf1/RNA polymerase II complex component, homolog (S. cerevisiae); Tsp; Tsbp; Sh2bp1; AA409336; mKIAA0155; RNA polymerase-associated protein CTR9 homolog; p150TSP; SH2 domain-binding protein 1; TPR-containing, SH2-binding phosphoprotein of 150 kDa; SH2 domain binding protein 1 (tetratricopeptide repeat containing); tetratricopeptide repeat-containing, SH2-binding phosphoprotein of 150 kDa;
Entrez Gene ID	22083
Protein Refseq	NP_033457
UniProt ID	Q05CJ7
Chromosome Location	7 F1; 7 57.98 cM
Pathway	PluriNetWork;
Function	SH2 domain binding; protein binding;