



Anti-LNX2 (aa 1-100) polyclonal antibody (DPAB-DC973)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	LNX2 (ligand of numb-protein X 2) is a protein-coding gene. Diseases associated with LNX2 include q fever, and endocarditis. GO annotations related to this gene include PDZ domain binding and zinc ion binding. An important paralog of this gene is PARD3B.
Immunogen	LNX2 (NP_699202, 1 a.a. ~ 100 a.a) partial recombinant protein with GST tag. The sequence is MGTTSDDEMVSVEQTSSSSLNPLCFECGQQHWTRENHLYNYQNEVDDDLVCHICLQPLLQP LDTPCGHTFCYKCLRNFLQEKDFCPLDRKRLHFKLCKKSS
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Cell lysate), WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name [LNX2 ligand of numb-protein X 2 \[Homo sapiens \(human\) \]](#)

Official Symbol	LNK2
Synonyms	LNK2; ligand of numb-protein X 2; PDZRN1; ligand of Numb protein X 2; numb-binding protein 2; PDZ domain containing ring finger 1; PDZ domain-containing RING finger protein 1;
Entrez Gene ID	222484
Protein Refseq	NP_699202
UniProt ID	Q8N448
Chromosome Location	13q12.2
Function	PDZ domain binding; protein binding; zinc ion binding;