



Anti-ZNRF2 (aa 146-207) polyclonal antibody (DPAB-DC976)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	ZNRF2 (zinc and ring finger 2) is a protein-coding gene. Diseases associated with ZNRF2 include neuronitis. GO annotations related to this gene include zinc ion binding and ligase activity. An important paralog of this gene is ZNRF1.
Immunogen	ZNRF2 (NP_667339, 146 a.a. ~ 207 a.a) partial recombinant protein with GST tag. The sequence is LPAHLSPHMFGGFKCPVCSKFVSSDEMDLHLMCLTKPRITYNEDVLSKDAGECAICLEE LQ
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name [ZNRF2 zinc and ring finger 2 \[Homo sapiens \(human\) \]](#)

Official Symbol	ZNRF2
Synonyms	ZNRF2; zinc and ring finger 2; RNF202; E3 ubiquitin-protein ligase ZNRF2; protein Ells2; RING finger protein 202; zinc finger/RING finger 2; zinc/RING finger protein 2;
Entrez Gene ID	223082
Protein Refseq	NP_667339
UniProt ID	A0A090N8Y8
Chromosome Location	7p14.3
Function	ligase activity; ubiquitin-protein transferase activity; zinc ion binding;