



Rabbit Anti-Pig CCL5 Polyclonal Antibody (DPAB-JX23111324)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	<p>Reacts with Recombinant Swine, Bovine, Ovine, Equine CCL5.</p> <p>Does not react with recombinant canine, chicken, cynomolgus monkey, dolphin, feline, ferret, guinea pig, hamster, human, mouse, turkey CCL5.</p>
Target	Swine CCL5
Immunogen	Recombinant Swine CCL5
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Swine
Purification	Affinity chromatography
Conjugate	Unconjugated
Applications	<p>ELISA, WB</p> <p>We recommend the following as antibody pair (Capture - Detection): DPAB-JX23111324 - DPAB-JX23111325</p> <p>Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.</p>
Format	Purified, Liquid
Concentration	Lot specific
Size	100 µg

Buffer	PBS
Preservative	None
Storage	Store at -20°C or below. Avoid multiple freeze-thaw cycles.

BACKGROUND

Introduction	CCL5, also known as RANTES (regulated and normal T cell expressed and secreted) is a member of the C-C Chemokine Family. There are at least 27 distinct members of the C-C subgroup reported for mammals. RANTES (CCL5) is chemotactic for T cells, eosinophils, and basophils, and plays an active role in recruiting leukocytes into inflammatory sites. CCL5 plays a role in many disease states, including rheumatoid arthritis, HIV, and cancer.
---------------------	---