



# Rabbit Anti-Zebrafish BAX Polyclonal Antibody (DPAB-JX2361606)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Target</b>	BAX
<b>Immunogen</b>	Peptide corresponding to Zebrafish BAX
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Zebrafish
<b>Purification</b>	Antigen Affinity Purification
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB: 1/500-1/1000 Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
<b>Format</b>	Purified, Liquid
<b>Concentration</b>	Lot specific
<b>Size</b>	50 µL, 100 µL
<b>Buffer</b>	PBS containing 50% glycerol, 0.5% BSA, pH 7.3.
<b>Preservative</b>	None
<b>Storage</b>	Store at 2-8°C for up to 2 weeks. Store at -20°C for long term storage in small aliquots. Avoid freeze/thaw cycles.

Ship

Wet ice

---

## BACKGROUND

### Introduction

Apoptosis regulator BAX, also known as bcl-2-like protein 4, is a protein that in humans is encoded by the BAX gene. BAX is a member of the Bcl-2 gene family.

---

## GENE INFORMATION

Entrez Gene ID

[58081](#)

---

UniProt ID

[Q9I9N4](#)

---