



# Rabbit Anti-Zebrafish CASP1 Polyclonal Antibody (DPAB-JX2361609)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Target</b>	Caspase 1
<b>Immunogen</b>	Peptide corresponding to Zebrafish caspase-1
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Zebrafish
<b>Purification</b>	Antigen Affinity Purification
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB: 1/500-1/1000 Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
<b>Format</b>	Purified, Liquid
<b>Concentration</b>	Lot specific
<b>Size</b>	50 µL, 100 µL
<b>Buffer</b>	PBS containing 50% glycerol, 0.5% BSA, pH 7.3.
<b>Preservative</b>	None
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Ship

Wet ice

---

## BACKGROUND

### Introduction

CASP1(caspase-1) is also named as IL1BC, IL1BCE and belongs to the peptidase C14A family.

---

## GENE INFORMATION

Entrez Gene ID

[57933](#)

---

UniProt ID

[Q9I9L7](#)

---