



Rabbit Anti-Zebrafish TLR3 Polyclonal Antibody (DPAB-JX2361661)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	TLR3
Immunogen	Peptide from zebrafish TLR3 sequence
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Zebrafish
Purification	Antigen Affinity Purification
Conjugate	Unconjugated
Applications	WB: 1/500-1/1000 Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Purified, Liquid
Concentration	Lot specific
Size	50 µL, 100 µL
Buffer	PBS containing 50% glycerol, 0.5% BSA, pH 7.3.
Preservative	None
Storage	Store at 2-8°C for up to 2 weeks. Store at -20°C for long term storage in small aliquots. Avoid freeze/thaw cycles.

Ship

Wet ice

BACKGROUND

Introduction

TLR3 (Toll-like receptor 3,CD283) is a member of the Toll-like receptor (TLR) family which plays a fundamental role in pathogen recognition and activation of innate immunity.

GENE INFORMATION

Entrez Gene ID

[403126](#)

UniProt ID

[Q32PW5](#)
