



Sheep Anti-Agkistrodon piscivorus Venom Polyclonal Antibody (DPAB-JX2371313)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Agkistrodon piscivorus Venom
Isotype	IgG
Source/Host	Sheep
Species Reactivity	N/A
Purification	Unpurified
Conjugate	Unconjugated
Applications	IA Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Unpurified, antiserum
Size	1 mL
Preservative	None
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Ship	Wet ice

45-1 Ramsey Road, Shirley, NY 11967, USA

Tel: 1-631-624-4882 Fax: 1-631-938-8221