

Sheep Anti-Vipera ammodytes Venom Polyclonal Antibody (DPAB-JX2371321)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| Immunogen | Vipera ammodytes Venom |
|--------------------|---|
| Isotype | IgG |
| Source/Host | Sheep |
| Species Reactivity | N/A |
| Purification | Unpurified |
| Conjugate | Unconjugated |
| Applications | IA Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded. |
| Format | Unpurified, antiserum |
| Size | 1 mL |
| Preservative | None |
| Storage | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |
| Ship | Wet ice |

© Creative Diagnostics All Rights Reserved