



## Sheep Anti-Notechis Venom Polyclonal Antibody (DPAB-JX2371336)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	Notechis Venom
Isotype	IgG
Source/Host	Sheep
Species Reactivity	N/A
Purification	Unpurified
Conjugate	Unconjugated
Applications	IA  Each laboratory should determine an optimum working titer for use in its particular application.  Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Unpurified, antiserum
Size	1 mL
Preservative	None
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Ship	Wet ice

Email: info@creative-diagnostics.com

Tel: 1-631-624-4882 Fax: 1-631-938-8221