



Sheep Anti-Human Calcitonin Polyclonal Antibody (DPAB-JXL23213)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|---------------------------|---|
| Specificity | There was no cross-reactivity with human katacalcin or human calcitonin gene related peptide |
| Immunogen | Synthetic human Calcitonin conjugated with BSA |
| Isotype | IgG |
| Source/Host | Sheep |
| Species Reactivity | Human |
| Purification | Unpurified |
| Conjugate | Unconjugated |
| Applications | RIA 1:80,000 Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded. |
| Reconstitution | Resuspend in 20 µl/100 µl aqua bidest. |
| Format | Unpurified, lyophilized |
| Size | 20 µl, 100 µl |
| Preservative | None |
| Storage | Store at 2-8 °C (lyophilized), -20 °C (dissolved) |

GENE INFORMATION

Synonyms

Calca; calcitonin-related polypeptide alpha; CAL6; CGRP; Cal1; Calc; RATCAL6; calcitonin; calcitonin gene-related peptide 1; CGRP-I; alpha-type CGRP; calcitonin gene-related peptide 1; calcitonin gene-related peptide I; calcitonin/calcitonin-related polypeptide, alpha

Entrez Gene ID

[796](#)

UniProt ID

[P01258](#)
