

Rabbit Anti-Human GLP-2 polyclonal antibody (DPAB-JXL23247)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the N-term region of human GLP-2
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Rat
Purification	Affinity chromatography
Conjugate	Unconjugated
Applications	ELISA 1:500-1:1,000; WB 1:200-1:500; IHC 1:200-1:500 Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Purified, Liquid
Concentration	1.0 mg/mL
Size	100 μg
Buffer	PBS pH7.4, 25% glycerol, 0.1 mg/ml BSA with 0.09% sodium azide
Preservative	0.09% sodium azide
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Ship	Wet ice

45-1 Ramsey Road, Shirley, NY 11967, USA

Email: info@creative-diagnostics.com

Tel: 1-631-624-4882 Fax: 1-631-938-8221

© Creative Diagnostics All Rights Reserved

BACKGROUND

Introduction

Glucagon-like peptide-2 (GLP-2) is a 33-aa peptide corresponding to aa 146-178 of glucagon. GLP-2 is liberated in the intestinal L cells after glucagon processing. GLP-2 stimulates intestinal growth and up-regulates villus height in the small intestine, concomitant with increased crypt cell proliferation and decreased enterocyte apoptosis. The gastrointestinal tract from the stomach to the colon is the principal target for the GLP-2 action. GLP-2 plays a key role in nutrient homeostasis, enhancing nutrient assimilation through enhanced gastrointestinal function as well as increasing nutrient disposal. GLP-2 also stimulates intestinal glucose transport and decreases mucosal permeability.

GENE INFORMATION

Synonyms	Glucagon; Glucagon-like peptide 2; GLP-2; GCG; GLP2
Entrez Gene ID	<u>2641</u>
UniProt ID	<u>P01275</u>