



Anti-HBV Surface Antigen Polyclonal antibody (DPAB-TJ021)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Goat Anti HBsAg Goat Antibody to Human Hepatitis B Surface Antigen (HBsAg)
Antigen Description	HBsAg is the surface antigen of the hepatitis B virus (HBV). It indicates current hepatitis B infection.
Specificity	Reacts with both ad and ay antigens. No reaction when tested against normal human serum and single arc when tested against HBsAg elevated human serum by IEP
Target	HBV Surface Antigen
Immunogen	Purified Human Hepatitis B Surface Ag.
Isotype	IgG
Source/Host	Goat
Species Reactivity	HBV
Conjugate	Unconjugated
Applications	ELISA, CLIA, RIA, IHC, WB
Format	Purified, Liquid
Concentration	2 mg/ml (OD280nm, E0.1% = 1.4)
Size	1 mg
Buffer	10MmPBS, pH 7.4

Preservative None

Storage Short term, store at 2-8°C. Long term, store at -20°C.

BACKGROUND

Introduction HBsAg is the surface antigen of the hepatitis B virus (HBV). It indicates current hepatitis B infection.

Keywords Anti-HBV Mab; HBsAg; Hepatitis B Virus Surface Antigen; Mouse Anti-Hepatitis B Surface Antigen Monoclonal Antibody; HBV; Hepatitis B Virus; HBV; Hepatitis B Virus; HBsAg; HBV major surface antigen; HBV surface antigen; Hepatitis B Virus major surface anti
