



Rabbit Anti-HBV Surface Antigen Polyclonal antibody [HRP] (DPAB0227-P)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Hepatitis B surface antigen purified from human serum. Mixture of subtypes ad & ay.
Isotype	IgG
Source/Host	Rabbit
Purification	Affinity purification
Conjugate	HRP
Applications	ELISA We recommend the following for sandwich ELISA (Capture - Detection): DMAB3526 - DPAB0227R DMAB3527 - DPAB0227R
Format	Liquid
Buffer	Phosphate Buffered Saline containing 10 mg/mL BSA, 0.002% Thimerosal
Storage	Store at -20°C.
Ship	Wet ice