



Anti- PLG polyclonal antibody (DPAB1475RC)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Chicken Plasminogen and plasmin
Immunogen	Chicken Plasminogen
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Purification	98% pure (SDS-PAGE). Protein A chromatography
Conjugate	Unconjugated
Applications	Suitable for use in ELISA (1:10,000) and Western blot (1-5ug/ml). Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Purified, Liquid
Concentration	8.5mg/ml (OD280nm, E0.1% = 1.36)
Size	1 mg
Buffer	0.05M Sodium phosphate, pH 6.6 containing 0.1M Sodium chloride and 1mM EDTA
Preservative	None
Storage	Store at -70°C. Avoid multiple freeze/thaw cycles.

GENE INFORMATION

Gene Name	PLG plasminogen [Gallus gallus]
Official Symbol	PLG
Synonyms	PLG; plasminogen; EC 3.4.21.7; DKFZp779M0222; plasmin; EC 3.4.21
Entrez Gene ID	421580
Protein Refseq	XP_419618
UniProt ID	F1NWX6
Chromosome Location	3
Pathway	Blood Clotting Cascade; Influenza A; Neuroactive ligand-receptor interaction