



Anti-native (native) polyclonal antibody (DPAB15093)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-native Polyclonal Antibody
Specificity	Specifically recognizes the light neurofilament subunit (~68 KDa).
Target	native
Immunogen	Native purified Nefl from pig spinal cord.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Birds, Bovine, Human, Mammals, Mouse, Pig, Rat
Conjugate	Unconjugated
Applications	WB, IHC, ICC, IF
Molecular Weight	68
Format	Liquid
Buffer	In antiserum
Preservative	None
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.