



Anti-ORM1 polyclonal antibody (DPAB2064RR)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit polyclonal to rat α -1-acid glycoprotein, affinity purified.
Antigen Description	An acute phase reactant protein; may have a role in the acute inflammation response.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Rat
Conjugate	Unconjugated
Applications	ELISA, WB
Concentration	0.5-2.5 mg/ml (Assuming A280 = 14 for a 10 mg/ml solution).
Size	100 μ g
Buffer	10 mM sodium phosphate, 150 mM NaCl, pH 7.2, 0.1 % NaN ₃ .
Preservative	0.1% Sodium Azide
Storage	2-8°C.

GENE INFORMATION

Gene Name	Orm1 orosomuroid 1 [Rattus norvegicus]
Official Symbol	Orm1
Synonyms	Orm1; orosomuroid 1; AGP; OMD; Orm; Agpa1; alpha-1-acid glycoprotein; ORM; AGP1; OTTHUMP00000022741; AGP 1; alpha-1-acid glycoprotein 1; OMD 1; orosomuroid-1; AGP-A;

Orosomucoid-1

Entrez Gene ID	24614
Protein Refseq	NP_445740
UniProt ID	P02764
Chromosome Location	5q24
Function	drug binding