



Anti-Salusin-beta polyclonal antibody (DPAB21488H)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|---------------------------|--|
| Product Overview | Rabbit Anti-Salusin-beta Polyclonal Antibody |
| Target | Salusin-beta |
| Immunogen | Salusin-beta (Human)-TGA |
| Isotype | IgG |
| Source/Host | Rabbit |
| Species Reactivity | Human |
| Conjugate | Unconjugated |
| Size | 50 µl |
| Preservative | None |
| Storage | -20 °C, Avoid freeze / thaw cycles |

BACKGROUND

| | |
|---------------------|---|
| Introduction | Salusin-beta - Undiluted Antiserum for Immunohistochemistry |
| Keywords | Salusin; Rabbit Anti-Salusin-beta Polyclonal Antibody; Anti-Salusin-beta Polyclonal Antibody; Salusin-beta Polyclonal Antibody Rabbit Anti-Salusin-beta PAb; Anti-Salusin-beta PAb; Salusin-beta PAb; Rabbit Anti-Salusin-beta Antibody; Anti-Salusin-beta Anti |