



Anti-Natrium+ Channel II polyclonal antibody (DPAB21637H)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Natrium+ Channel II Polyclonal Antibody
Specificity	Humanin 100%, Humanin (5-24) 35%, Humanin (9-24) 1.9%, Humanin (14-24) 0.18%, Humanin (19-24) 0.18%, Humanin (9-18) 0%, Rattin-38 0%, Rattin-24 0%.
Target	Natrium+ Channel II
Immunogen	Synthetic Natrium+ Channel II (Rat, 467-485) carrier free
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Rat
Conjugate	Unconjugated
Applications	IHC
Size	50 µl
Preservative	None
Storage	-20 °C, Avoid freeze / thaw cycles