



Anti-Rattin-24 polyclonal antibody (DPAB21641H)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Rattin-24 Polyclonal Antibody
Specificity	Rattin-24 100% Rattin-38 100% Humanin 3.8%
Target	Rattin-24
Immunogen	Synthetic Rattin-24-KLH
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Rat
Conjugate	Unconjugated
Applications	RIA, IHC
Size	50 µl
Preservative	None
Storage	-20 °C, Avoid freeze / thaw cycles

BACKGROUND

Introduction	This antiserum was raised in a rabbit by immunization with a carrier free synthetic Rattin-24 peptide. The product vial contains 50 µL of the titled antiserum obtained by lyophilizing its 0.001 M phosphate buffer (pH 7.0, 0.5mL) solution. It can be used for immunoassay, immunohistochemistry, and tec.
--------------	---

Keywords

Rattin-24 antibody; Rabbit Anti-Rattin-24 Polyclonal Antibody; Anti-Rattin-24 Polyclonal Antibody; Rattin-24 Polyclonal Antibody Rabbit Anti-Rattin-24 PAb; Anti-Rattin-24 PAb; Rattin-24 PAb; Rabbit Anti-Rattin-24 Antibody; Anti-Rattin-24 Antibody; Rattin-
