



# Rabbit Anti-Streptococcus Group A Polyclonal antibody (DPAB2782)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Product Overview	Rabbit polyclonal to bacteria Group A Streptococcus.
Target	Streptococcus Group A
Immunogen	Bacteria, Variable
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Streptococcus group A
Conjugate	Unconjugated
Applications	IP, ELISA, WB, IF
Size	0.5 ml
Preservative	0.1% Sodium Azide
Storage	Store at 2-8 degrees C, For long term storage store frozen.

## BACKGROUND

Introduction	The genus Streptococcus is comprised of a wide variety of both pathogenic and commensal gram positive bacteria which are found to inhabit a wide range of hosts, including humans, horses, pigs and cows. Within the host, streptococci are often found to colonize the mucosal surfaces of the mouth, nose and pharynx. Streptococci can be divided into many groups on the basis of antigenic differences in group-specific polysaccharides located in the bacterial cell wall.
--------------	---

**Keywords**

GAS; Streptococcus pyogenes

---