



Anti-Staphylococcus Group A Polyclonal antibody (DPAB2783)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Goat polyclonal to bacteria Group A Streptococcus.
Target	Staphylococcus Group A
Immunogen	Bacteria, Variable
Source/Host	Goat
Species Reactivity	Staphylococcus
Conjugate	Unconjugated
Applications	IP, ELISA, WB, IF
Size	0.5 ml
Preservative	None
Storage	Store at 2-8 degrees C, For long term storage store frozen.

BACKGROUND

Introduction	Streptococcus pyogenes is a spherical, Gram-positive bacterium that is the cause of group A streptococcal infections. S. pyogenes displays streptococcal group A antigen on its cell wall. S. pyogenes typically produces large zones of beta-hemolysis (the co
Keywords	GAS; Streptococcus pyogenes