



Anti-Gamma Globulins polyclonal antibody (DPAB2792)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit polyclonal to guinea pig Gamma Globulins.
Immunogen	Guinea Pig, Variable
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Guinea pig
Conjugate	Unconjugated
Applications	IP, ELISA, WB, IF
Size	2 ml
Preservative	None
Storage	Store at 2-8 degrees C, For long term storage store frozen.

BACKGROUND

Introduction	Gamma globulins are a class of globulins, identified by their position after serum protein electrophoresis. The most significant gamma globulins are immunoglobulins ("Igs"), more commonly known as antibodies, although some Igs are not gamma glob
Keywords	IFNG; 310C; AMCF1; Beta thromboglobulin like protein; CXC chemokine ligand 8; CXCL8; Emotakin; GCP1; Granulocyte chemotactic protein 1; IFG; IFI; IFN Gamma; IFN Immune; IFN-gamma; IFNG; IFNG_HUMAN; Immune Interferon; Inteleukin 8; Interferon gamma; Inter