



Anti-Ganglioside GM2 polyclonal antibody (DPAB2897)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Ganglioside GM2 Polyclonal Antibody, Recognizes ganglioside GM2.
Specificity	a broad range of species
Immunogen	purified, bovine ganglioside GM2
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Bovine
Conjugate	Unconjugated
Applications	ELISA, IB
Format	Liquid. Undiluted serum.
Size	100 µl
Preservative	None
Storage	-20°C, Avoid freeze/thaw
Ship	Blue Ice Only Standard Handling

BACKGROUND

Introduction	GM2 is a type of ganglioside. It is associated with GM2 gangliosidoses such as Tay-Sachs disease.
--------------	---------------------------------------------------------------------------------------------------

Keywords

Ganglioside GM2; GM2 (ganglioside); GM2; Ganglioside GM2; Tay-Sachs ganglioside; GANGLIOSIDE GM2 RABBIT; NGNA; β -D-GalNAc-(1 \rightarrow 4)-[α -Neu5Ac-(2 \rightarrow 3)]- β -D-Gal-(1 \rightarrow 4)- β -D-Glc-(1 \leftrightarrow 1)-N-octadecanoylsphingosine; GANGLIOSIDE GM2
