



# Anti-Pneumadin polyclonal antibody (DPAB3226)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Polyclonal Antibody to Pneumadin
<b>Specificity</b>	Synthetic human Pneumadin
<b>Immunogen</b>	Synthetic human Pneumadin (AGEPKLDAGV), poly Lysin conjugated
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Conjugate</b>	Unconjugated
<b>Format</b>	serum
<b>Concentration</b>	20 µl / 100 µl (lyophilized) resuspend in 20 µl / 100 µl aqua bidest
<b>Preservative</b>	None
<b>Storage</b>	2°C-8°C (lyophilized); - 20°C (dissolved) Repeated thawing and freezing must be avoided

## BACKGROUND

<b>Introduction</b>	Pneumadin, a new lung peptide which triggers antidiuresis is a rabbit polyclonal antibody to protein/bioactive peptide. Small volumes of Pneumadin antibody may occasionally become entrapped in the seal of the product vial during shipment and storage. If necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the container's cap. Certain products may require to ship with dry ice.
---------------------	---

**Keywords**

Pneumadin; PNEU; PNM

---