



Anti-Benzoylecgonine polyclonal antibody (DPAB3466)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Sheep Polyclonal antibody to Benzoylecgonine (Cocaine).
Immunogen	Benzoylecgonine - Thyroglobulin
Source/Host	Sheep
Species Reactivity	N/A
Conjugate	Unconjugated
Format	Neat serum
Preservative	0.1% Sodium Azide
Storage	Store at 2-8°C

BACKGROUND

Introduction	Benzoylecgonine is the predominant cocaine metabolite. This compound is not biologically active but is used as an indicator of cocaine use.
Keywords	Benzoylecgonine; Cocaine metabolite; BENZOYLECGONINE; O-Benzoylecgonin; 1alphaH,5alphaH-Tropane-2beta-carboxylic acid, 3beta-hydroxy-, benzoate (ester); 3-(Benzoyloxy)-8-methyl-8-azabicyclo[3.2.1]octane-2-carboxylic acid; 8-Azabicyclo[3.2.1]octane-2-carbo