



Anti-Gentamycin polyclonal antibody (DPAB3473)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Polyclonal antibody to Gentamycin antigen.
Immunogen	Gentamycin
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Conjugate	Unconjugated
Format	Neat serum.
Size	250 µl
Preservative	None
Storage	2-8°C.

BACKGROUND

Introduction	Gentamicin is an aminoglycoside antibiotic, used to treat many types of bacterial infections, particularly those caused by Gram-negative organisms. However, gentamicin is not used for <i>Neisseria gonorrhoeae</i> , <i>Neisseria meningitidis</i> or <i>Legionella pneumophila</i> . Gentamicin is also ototoxic and nephrotoxic, with this toxicity remaining a major problem in clinical use.
Keywords	Gentamicin; GENTAMICIN; GENTAMYCIN; GENTAMYCINE; GentamysinsolutionforBiochemistry; Gentacycol; Gentavet; GENTAMICINUM; Refobacin tm;

Uromycine; Gentamicin (base and/or unspecified salts); Lyracylin; 2-[4,6-Diamino-3-[3-amino-6-(1-methylaminoethyl)oxan-2-
