



Anti-Triiodothyronine polyclonal antibody (DPAB3484)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Goat Polyclonal antibody to human Triiodothyronine antigen.
Target	Triiodothyronine
Immunogen	Triiodothyronine-BSA
Source/Host	Goat
Species Reactivity	Human
Conjugate	Unconjugated
Format	Neat serum
Size	1 ml
Preservative	None
Storage	2-8°C short term; -20°C long term.

BACKGROUND

Introduction	Triiodothyronine, C ₁₅ H ₁₂ I ₃ NO ₄ , also known as T ₃ , is a thyroid hormone. It affects almost every physiological process in the body, including growth and development, metabolism, body temperature, and heart rate.
Keywords	Triiodothyronine; Liothyronine; T ₃ ; (2S)-2-amino-3- [4-(4-hydroxy-3-iodo-phenoxy)- 3,5-diiodo-phenyl]propanoic acid; LIOETHYRONIN; L-3,3",5"-TRIIODOTHYRONINE; (3,5,3")-TRIIODO-L-THYRONINE; Reverse T ₃ ; 3,3",5"-triiodo-L-thyronine free acid; Isoliothyronine