



Anti-Testosterone polyclonal antibody (DPAB3485)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Polyclonal antibody to human Testosterone antigen.
Immunogen	Testosterone-19
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Conjugate	Unconjugated
Format	Neat serum
Preservative	None
Storage	4-8°C short term; -20°C long term.

BACKGROUND

Introduction	Testosterone is a steroid hormone from the androgen group and is found in mammals, reptiles, birds, and other vertebrates. In mammals, testosterone is primarily secreted in the testes of males and the ovaries of females, although small amounts are also secreted by the adrenal glands. It is the principal male sex hormone and an anabolic steroid.
Keywords	Testosterone; 17b hydroxy 4 androsten 3 one; testred; TESTOSTERONE; Testex; testoderm; TESTOSTERON; ORQUISTERON; oreton; PRIMOTESTON; virilon; component of Di-Genik; component of Di-Met; component of Duogen; component of Mal-O-Fem; Cristerona T; cristeron