



Rabbit Anti-H12N5 (A/pintail/Alaska/102/05) HA Polyclonal Antibody (DPAB3961)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Polyclonal antibody to H12N5 Hemagglutinin.
Specificity	Reacts with hemagglutinin (H12N5) protein. Cross reactivity to other subtypes not tested.
Target	IAV H12N5 Hemagglutinin
Immunogen	in vivo expressed hemagglutinin (H12N5)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	IAV
Purification	Immunoaffinity chromatography
Conjugate	Unconjugated
Applications	WB, ELISA
Size	100 µg
Preservative	None
Storage	Store at 4°C; DO NOT FREEZE; Stable for 6 months from the date of shipment. Non-hazardous.

BACKGROUND

Introduction	H12N5 is a subtype of the species Influenza A virus (bird flu virus)
---------------------	--

Keywords

HA; HA1; HA2; Hemagglutinin; Hemagglutinin HA1 chain; Hemagglutinin HA2 chain; H12N5 hemagglutinin; H12 (H12N5); Group V; Orthomyxoviridae; Influenzavirus A; Influenzavirus B; Influenzavirus C; Isavirus; Thogotovirus
