



Rabbit Anti-H15N8(A/Australia/83/2000) HA Polyclonal Antibody (DPAB3964)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Polyclonal antibody to H15N8 H15.
Specificity	Reacts with the H15 protein of influenza A virus. Cross-reactivity to other subtypes not tested.
Target	IAV H15N8 H15
Immunogen	H15 (H15N8)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	IAV
Purification	Immunoaffinity chromatography
Conjugate	Unconjugated
Applications	WB
Size	100 µg
Preservative	None
Storage	Store at 4°C; DO NOT FREEZE; Stable for 6 months from the date of shipment. Non-hazardous.

BACKGROUND

Introduction	H15 is a subtype of the species Influenza A virus (bird flu virus).
--------------	---

Keywords

H15 (H15N8); H15N8 H15; Group V; Orthomyxoviridae; Influenzavirus A; Influenzavirus B; Influenzavirus C; Isavirus; Thogotovirus
