



Rabbit Anti-H7N7 (A/Chicken/Netherlands/1/03) NA Polyclonal Antibody (DPAB3966)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Polyclonal antibody to H7N7 NA.
Specificity	Reacts with Influenza A NA (H7N7) protein. Cross-reactivity with other NA proteins not tested.
Target	IAV H7N7 Neuraminidase
Immunogen	in vivo expressed NA from Influenza A virus (H7N7)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	IAV
Purification	Immunoaffinity chromatography
Conjugate	Unconjugated
Applications	WB, ELISA
Size	100 µg
Preservative	None
Storage	Store at 4°C; DO NOT FREEZE; Stable for 6 months from the date of shipment. Non-hazardous.

BACKGROUND

Introduction	H7N7 is a subtype of the species Influenza A virus (bird flu virus).
---------------------	--

Keywords

H7N7 NA; NA (H7N7); Group V; Orthomyxoviridae; Influenzavirus A; Influenzavirus B; Influenzavirus C; Isavirus; Thogotovirus
