



Goat anti-Glutathione-S-Transferase (GST) Polyclonal antibody (DPAB4081)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Glutathione-S-Transferase (GST) Affinity Purified Antibody, HRP Conjugated
Immunogen	GST from Schistosoma japonicum coupled to KLH.
Source/Host	Goat
Conjugate	Unconjugated
Applications	ELISA, WB, IHC/ICC
Size	1 ml
Buffer	50mM HEPES, pH 7.1, 0.1M NaCl
Preservative	None
Storage	Store at 4°C; stable for one year

BACKGROUND

Introduction	GST (Glutathione S-Transferase) is a 26kDa protein encoded by the parasitic helminth Schistosoma japonicum and widely used in the pGEX family of GST plasmid expression vectors as a fusion protein with foreign proteins.
Keywords	Glutathione S Transferase; Glutathione S transferase Mu 1; GST 26; GST class mu 1; GST1; GSTM1 1; GSTM1a 1a; GSTM1b 1b; GTH4; GTM1; H B; HB subunit 4; MGC26563; MU 1; MU; SJ26 antigen

GENE INFORMATION

Entrez Gene ID [2944](#)

UniProt ID [P09488](#)
