



Anti-Sulfametopyrazine polyclonal antibody (DPAB4208)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit polyclonal antibody to sulfametopyrazine.
Immunogen	Sulfametopyrazine-KLH Conjugate
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Conjugate	Unconjugated
Buffer	0.01 mol/L pH 7.4 Phosphate Buffer Solution.
Preservative	None
Storage	Short-term store at 2-8°C, Long-term store at -20°C.

BACKGROUND

Introduction	Sulfalene is a sulfonamide antibacterial used for the treatment of chronic bronchitis, urinary tract infections and malaria.
Keywords	4-amino-n-(3-methoxypyrazinyl)benzenesulfonamide; as18908; dalysep; f.i.5978; farmitalia204/122; kelfizin; kelfizina; kelfizine; longum; n(sup1)-(3-methoxy-2-pyrazinyl)-sulfanilamid; n(sup1)-(3-methoxy-2-pyrazinyl)sulfanilamide; n(sup1)-(3-methoxypyrazinyl)-sulfanilamid; n(sup1)-(3-methoxypyrazinyl)sulfanilamide; nsc-110433; polycidal; smp2; sulfamethopyrazine; sulfametopyrazine; sulfapyrazinemethoxine; sulfapyrazinemethoxyine; sulfapyrazinemethoxyne