



Anti-Sulfadoxine polyclonal antibody (DPAB4214)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit polyclonal antibody to sulfadoxine.
Immunogen	Sulfadoxine-KLH Conjugate
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Conjugate	Unconjugated
Buffer	0.01 mol/L pH 7.4 Phosphate Buffer Solution.
Preservative	None
Storage	Short-term store at 2-8°C, Long-term store at -20°C.

BACKGROUND

Introduction	Sulfadoxine (also spelled sulphadoxine) is an ultra-long-lasting sulfonamide often used in combination with pyrimethamine to treat or prevent malaria. It is also used, usually in combination with other drugs, to treat or prevent various infections in livestock.
Keywords	4-amino-n-(5,6-dimethoxy-4-pyrimidinyl)-benzenesulfonamid; 6-(4-aminobenzenesulfonamido)-4,5-dimethoxypyrimidine; fanasil; fanzil; n'-(5,6-dimethoxy-4-pyrimidyl)sulfanilamide; ro4-4393; sulphadoxine; sulphormethoxine; wr4873; SULFADOXIN; SULFADOXINE; sulformethoxine; sulforthomidine; 6-(4-aminobenzenesulfonamide)-4,5-dimethoxypyrimidine; 4-sulfanilamido-5,6-dimethoxypyrimidine