



# Anti-Cathepsin L Protease polyclonal antibody (DPAB4239)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Rabbit polyclonal antibody to cathepsin L protease [proregion].
<b>Specificity</b>	Antiserum reacts specifically with 36 kDa wild-type protein of C.elegans.
<b>Immunogen</b>	Synthetic peptide derived from Pro-region of CPL-1.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	N/A
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	IEM, WB
<b>Format</b>	Antiserum
<b>Size</b>	100 µl
<b>Preservative</b>	None
<b>Storage</b>	Short term storage: +4°C. Long term storage: -20°C

## BACKGROUND

<b>Introduction</b>	Cathepsin L can refer to: Cathepsin L1, previously called cathepsin L. Cathepsin L2 or cathepsin V.
<b>Keywords</b>	Cathepsin L1 light chain; CathepsinL; CATL; CATL1_HUMAN; CTSL; CTSL1; Major excreted

