



Rabbit Anti-Zebrafish Casp3b polyclonal antibody (DPABB-JX318)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide derived from the C-terminal region of zebrafish Caspase-3b p12 domain
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Zebrafish
Conjugate	Unconjugated
Applications	WB 1:500-1:1000 Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Liquid
Concentration	Lot specific
Size	50 µl
Buffer	MOPS buffer, 0.1% BSA, 50% Glycerol
Preservative	0.05% Sodium azide
Storage	For short-term storage, store at 2-8°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

BACKGROUND

Introduction

Caspase-3 is a caspase protein that interacts with caspase-8 and caspase-9. It is encoded by the CASP3 gene. CASP3 orthologs have been identified in numerous mammals for which complete genome data are available.

Keywords

Caspase 3, Apoptosis-Related Cysteine Peptidase B , Casp3L , Casp3B
