



Rabbit Anti-Human TRPV1 polyclonal antibody (DPABC-JX155)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	Vanilloid R1/TRPV1
Immunogen	Synthetic peptide corresponding to human VR1
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ICC/IF, IHC Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Liquid
Concentration	Lot specific
Size	100 µg
Buffer	PBS with 1 mg/ml BSA
Preservative	0.05% Sodium Azide
Storage	Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

GENE INFORMATION

Entrez Gene ID [7442](#)

UniProt ID [Q8NER1](#)
